



No hydration required

Hydrofera Blue Ready foam

Stock No	Product Size	Quantity/Box	HCPCS*
Standard Dressings			
HBRS2520	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRS4520	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRS8820	8" x 8" (20.3 cm x 20.3 cm)	10	A6211



Hydrate before use

Hydrofera Blue foam

Stock No	Product Size	Quantity/Box	HCPCS*
Standard Dressings			
HB2214	2" x 2" (5 cm x 5 cm)	10	A6209
HB4414	4" x 4" (10 cm x 10 cm)	10	A6209
HB6614	6" x 6" (15 cm x 15 cm)	10	A6210
HBT0906	9-mm Tunneling Dressing (1.2 g)	10 pair	A6215
Heavy Drainage Dressings			
HBHD4450	4" x 4" x .5" thick (10 cm x 10 cm x 1.2 cm)	5	A6209
HBHD6675	6" x 6" x .75" thick (15 cm x 15 cm x 1.9 cm)	5	A6210
Dressings With Moisture-Retentive Film			
HBRF4414	4" x 4" (10 cm x 10 cm)	10	A6209
HBRF2258	2.25" x 8" (5.7 cm x 20 cm)	10	A6210
HBRF4475	4" x 4.75" (10 cm x 12 cm) Island Dressing; <i>pad size 2" x 2.75" (5 cm x 6.5 cm)</i>	10	A6212
HBRF2650	2.5" (6 cm) Diameter Ostomy Dressing	10	A6209

1. Data on file, EN 13726-1; 2002.
2. Based on a January 2014 analysis of publicly available pricing of major brand non-adhesive, silver foams available on the Internet as compared to the per-unit suggested retail price of Hydrofera Blue Ready.

Manufactured for
Hollister Incorporated
2000 Hollister Drive
Libertyville, IL 60048 USA

Distributed in Canada by:
Distribué au Canada par:
Hollister Limited
95 Mary Street
Aurora, Ontario L4G 1G3

www.hollisterwoundcare.com

Rx Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

For product questions, sampling needs, or detailed clinical questions concerning our products in the US call **1.888.740.8999**. In Canada call **1.800.263.7400**.

*HCPCS are for reference only and subject to change.

The Hollister logo is a trademark of Hollister Incorporated. Hydrofera Blue Ready and the Hydrofera Blue logo are trademarks of Hydrofera, LLC.

All other trademarks are the property of their respective owners.

Made and printed in the USA.

©2014 Hollister Incorporated.

922496-0214

Hydrofera Blue Ready. Set. Go.



Introducing **Hydrofera Blue Ready** foam

- Antibacterial properties to help prevent infection
- No hydration required prior to application
- No secondary dressing required
- Up to 7-day wear time

When managing bioburden is a priority,
Hydrofera Blue Ready foam is the convenient choice.

Waterproof backing


- Provides a barrier to moisture and bacteria

Broad-spectrum antibacterial protection


- Provided by two organic pigments, methylene blue and gentian violet

Absorption

- Effectively absorbs exudate away from the wound



Hydrofera Blue Ready foam has absorption comparable to other antimicrobial foam dressings, based on *in vitro* testing¹



(Error bars represent one standard deviation)

Additional benefits

- Costs less than most silver foam dressings²
- Non-cytotoxic

Hydrofera Blue Ready foam.
 Ready. Set. Go.

Ready without hydration.
No secondary dressing required; just secure.
Good to go for up to 7 days.

Application

Ready

Select and apply the appropriate dressing size to ensure the dressing will cover the entire wound. With **Hydrofera Blue Ready** foam, no hydration is required.



Set

Secure the dressing using a bandage or tape.



Go

Dressing can be left in place up to 7 days. Can be used in conjunction with compression wraps and total contact casting.



Hydrofera Blue Ready foam has broad-spectrum antibacterial activity against microorganisms commonly found in wounds.

BACTERIA

- | | |
|---|---------------------------------|
| Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) | <i>Streptococcus pyogenes</i> |
| <i>Staphylococcus aureus</i> (coagulase negative) | <i>Klebsiella pneumoniae</i> |
| <i>Pseudomonas aeruginosa</i> | <i>Proteus mirabilis</i> |
| <i>Pseudomonas fluorescens</i> | <i>Proteus vulgaris</i> |
| <i>Escherichia coli</i> | <i>Enterobacter aerogenes</i> |
| <i>Enterococcus faecalis</i> | <i>Bacillus subtilis</i> |
| Vancomycin-resistant <i>Enterococcus</i> (VRE) | <i>Yersinia enterocolitica</i> |
| | <i>Serratia marcescens</i> |
| | <i>Streptococcus pneumoniae</i> |
| | <i>Enterobacter cloacae</i> |

YEAST

- Candida albicans*
- Candida krusei*
- Candida glabrata*

Indications for Use

Indicated for acute, traumatic, and chronic wounds, such as:

- | | |
|----------------------|------------------------|
| Pressure ulcers | Donor sites |
| Diabetic ulcers | Postsurgical incisions |
| Venous stasis ulcers | Trauma wounds |
| Arterial ulcers | Abrasions |
| Superficial burns | Lacerations |

